

# Microfabricated Polymeric Vessel Mimetics for 3D Cancer Cell Culture in Matrigel



Yiqi Cao<sup>1,2</sup>, Joanna E Poprawski<sup>2</sup>, Davi da Silva<sup>1</sup>, Matthew D Hall<sup>2</sup>, Randall H Pursley<sup>3</sup>, Marcial A Garmendia-Cedillos<sup>3</sup>, Nicole Y Morgan<sup>1</sup>, Thomas J Pohida<sup>3</sup> and Michael M Gottesman<sup>2</sup>

<sup>1</sup>National Institute of Biomedical Imaging and Bioengineering, <sup>2</sup>National Cancer Institute, <sup>3</sup>Center for Information Technology, National Institutes of Health, Bethesda, Maryland 20892, USA

## INTRODUCTION

Delivering a new drug to market in the United States costs an average of \$1 billion dollars over 15 years<sup>1</sup>. The failure rate of drugs entering early clinical trials is about 85%<sup>2</sup>. The divergence between results from preclinical studies and efficacy in human clinical trials is a persistent challenge in drug discovery.


**Animal models** are expensive, time consuming, and inadequate in predicting clinical efficacy  
**2D monolayer cell culture** fails to recapitulate the biological cues and cell-cell interaction found in the *in vivo* tumor microenvironment

**Standard 3D scaffolds and matrices** lack vasculature, thus limiting the transport of O<sub>2</sub> and other nutrients necessary for sustaining metabolism and growth.  
Here we report a 3D cell culture system using a microfabricated polymeric vessel mimetic that facilitates O<sub>2</sub> delivery to cancer cell cultures in anoxia. This study has implications for the development of an *in vitro* tumor models that mimic the tumor microenvironment.


## OBJECTIVES

- Aim 1: Bioreactor design.** Use microfabricated polydimethylsiloxane (PDMS) pillars to “vascularize” 3D cell cultures and mimic *in vivo* oxygen perfusion
- Aim 2: Validate cell culture growth pattern.** Evaluate the effect of the oxygen gradient on the growth patterns of tumor cell cultures
- Aim 3: Culture cells for drug screening.** Develop and optimize techniques to use high-throughput, quantitative cytotoxicity assays for 3D cultures.


## 1<sup>st</sup> GENERATION BIOREACTOR




Single chamber bioreactor system (11.5 cm x13 cm): two acrylic chambers separated by membrane with pillars (either PDMS or silicone hydrogel). Continuous gas flow maintains each chamber at its desired oxygen concentration<sup>3</sup>.



Micropillars for oxygen delivery are fabricated using template of SU-8 micropillars with 100 μm diameter and 250 μm height. Using the SU-8 master, negative PDMS molds were replicated, and silicone hydrogel micropillars were cast from PDMS molds. Dimensions and spacing of pillars are designed to mimic those of capillaries *in vivo*<sup>3</sup>.




OVCar8-dsRed2 cells with Hoechst stain after 7 d culture on PDMS (top) and silicone hydrogel (bottom) pillars. PDMS replaced SiHy as the pillar material in later design due to improved batch-to-batch consistency of physical properties in fabrication. Scale bars are 100 μm.




OVCar8-dsRed2 cells grown for 7 days on silicone hydrogel pillars with a 3% O<sub>2</sub> gradient were stained with hypoxia-marking dye pimonidazole and imaged (left). Cells become hypoxic about 100 μm from the pillar (circled in red), comparable to *in vivo* results<sup>3</sup>.

## 2<sup>nd</sup> GENERATION BIOREACTOR




The bioreactor is designed using SolidWorks and fabricated by laser cutting acrylic sheets. The top chamber and bottom chamber were sealed by gas-tight gaskets.

The redesigned bioreactor chamber accommodates 6 custom six-well cell culture plates for high-throughput experiments. The system features feedback control for the O<sub>2</sub> concentration in the hypoxic chamber and real-time monitoring of O<sub>2</sub> levels.




Custom six-well plate assembly: two acrylic plates in the six-well plate format secure the PDMS membrane in the middle using rare-earth magnets.

Vertical cross-section of a single well on the PDMS membrane. Each well contains an inset miniwell with microfabricated micropillars lining the bottom, where cells in Matrigel are seeded.




Vertical cross-section of a single miniwell containing micropillars. PDMS replaced silicone hydrogel as the material for the pillars because of increased batch-to-batch consistency in fabrication. Cells are embedded in Matrigel and seeded over the pillars, which deliver oxygen to the 3D cell layer.

## 3D CELL CULTURE IN BIOREACTOR




### The Effect of Plasma Cleaning PDMS Membranes



While plasma cleaning renders native PDMS hydrophilic and improved the homogenous deposition of Matrigel on the PDMS membrane, plasma cleaning also caused preferential cell adhesion on PDMS surfaces on the bottom of the membrane and on the pillars, thus confounding results. Plasma cleaning was therefore discontinued. Scale bars are 200 μm.

### Fluorescent Tracking of Cancer Cell Lines




CellTracker (Life Technologies), a fluorescent dye, was adapted for cell cultures in Matrigel and was used to visualize OvCar8 spheroids for live cell tracking. This technique can be used to visualize other non-fluorescent cell lines for culture in the bioreactor. Scale bar is 400 μm.

### Confocal Imaging of Custom Plates for Longitudinal Studies




A custom stage adaptor was designed in SolidWorks and 3D printed, to allow for confocal imaging of the intact assembly using a standard microscope. Keeping the stage intact enables longitudinal monitoring, instead of endpoint imaging only.

## CYTOTOXICITY ASSAY



	KB 3-1	CP.5		KB 3-1	CP.5		KB 3-1	CP.5
IC <sub>50</sub> (μM)	2.480	21.031	IC <sub>50</sub> (μM)	4.471	22.921	IC <sub>50</sub> (μM)	9.356	32.133
SD (μM)	0.102	2.831	SD (μM)	0.702	4.633	SD (μM)	0.702	6.868

A protocol was developed for performing cytotoxicity assays in cells cultured in Matrigel using CellTiter-Glo®. The assay shows that, as expected, cell resistance is escalated for multicellular spheroids as compared to single cells. Dose response curves of Kb3-1 cells, a parental cell lines, and Cp.5, a cisplatin-resistant cell line, are shown above for cisplatin, a platinum based anti-cancer drugs.



## FUTURE WORK

### I. Further characterization of the bioreactor system

- Growth pattern and cell localization
- RNA expression
- Drug sensitivity/penetration
- Hypoxia markers

### II. Culture other cancer cell lines in bioreactor

- Breast – MCF7
- Ovarian- OvCar 8
- Kidney-SN12C
- AML - CEM
- Liver- FOCUS
- Colon- DLD1

## ACKNOWLEDGMENTS

This research was funded by the Intramural Research Programs at the National Institutes of Biomedical Imaging and Bioengineering, the National Cancer Institute, and the Center for Information Technology at the National Institutes of Health.



## REFERENCES

- Kessel, M. The problems with today's pharmaceutical business--an outsider's view. *Nat Biotech* **29**, 27–33 (2011).
- Ledford, H. Translational research: 4 ways to fix the clinical trial. *Nature*, 477 (2011), pp. 526–528
- Jaeger AA, Das CK, Morgan NY, et al. Microfabricated polymeric vessel mimetics for 3-D cancer cell culture. *Biomaterials*. 2013;34(33):8301-13.